

Characterization of the Newly Designed Wall-Free Particle Evaporator (WALL-E) for Online Measurements of Atmospheric Particles

Linyu Gao^{1,†}, Imad Zgheib^{1,#,†}, Evangelos Stergiou², Cecilie Carstens¹, Félix Sari Doré^{1,‡},
Michel Dupanloup¹, Frederic Bourgain¹, Sébastien Perrier¹, Matthieu Riva^{1*}

¹ Univ Lyon, Université Claude Bernard Lyon 1, CNRS, IRCELYON, F-69626, Villeurbanne,
France

² Environmental Chemical Processes Laboratory (ECPL), Department of Chemistry, University
of Crete, Voutes Campus, 70013 Heraklion, Greece

[#] now at TOFWERK, 3645 Thun, Switzerland

[‡] now at Department of Chemistry and Molecular Biology, Atmospheric Science, University of
Gothenburg, SE-41390, Gothenburg, Sweden

[†] These authors contributed equally to this work

Corresponding to: matthieu.riva@ircelyon.univ-lyon1.fr

Abstract

Organic aerosols (OA) play a critical role in the atmosphere by directly altering human health and climate. Understanding their formation and evolution as well as their physicochemical properties requires a detailed characterization of their chemical composition. Despite advanced analytical techniques developed within the last decades, real-time online measurement of atmospheric particles remains challenging and suffers from different artifacts (i.e., thermal decomposition, fragmentation, wall loss). In this work, we introduce the newly designed wall-free particle evaporator (WALL-E) coupled with a chemical ionization mass spectrometer (CIMS) using bromide (Br^-) as the reagent ion. We comprehensively evaluate the performance of the WALL-E system, demonstrating its ability to evaporate particles while maintaining the integrity of the compounds composing the particles (i.e., minimal thermal decomposition). To demonstrate WALL-E's performance, the composition of aerosol particles formed from α -pinene ozonolysis in the presence of SO_2 is characterized. In addition, by applying the scan declustering method, we can now provide a quantification of the different species present in the condensed phase, e.g., $\text{C}_{10}\text{H}_{16}\text{O}_4$ 84 ng m^{-3} , $\text{C}_{19}\text{H}_{28}\text{O}_7$ 7 ng m^{-3} for a total SOA mass of 1 $\mu\text{g m}^{-3}$. While dimers exhibit higher sensitivities, they account for only 14-18% of the total particle masses, which is considerably lower than their signal fractions (23-29%). This suggests a potential overestimation of the dimer contributions when relying solely on signal fractions. In addition, volatility analysis using thermograms reveals a clear relationship between T_{50} and compound saturation vapor pressure (C^*), with lower-volatility species desorbing at higher temperatures. In addition, the measured T_{50} (the temperature at which 50% of a compound evaporates), for α -pinene-derived SOA products agree well with theoretical volatility estimation models (e.g., SIMPOL). Overall, this study demonstrates that WALL-E system coupled to a CIMS is a promising technique for real-time particle characterization (i.e., composition, quantification, and volatility) of atmospheric aerosols.

1 Introduction

Atmospheric organic aerosol (OA) particles play a critical role in the Earth's climate system and atmospheric processes by affecting the radiative forcing, cloud formation and albedo, atmospheric chemistry, environmental sustainability, and human health (Fehsenfeld et al., 1992; Laothawornkitkul et al., 2009; Mellouki et al., 2015; Charnawskas et al., 2017). A large fraction (20-90 %) of fine particles are comprised of OA (Kanakidou et al., 2005), which are estimated to have a global source of 150 Tg yr⁻¹ (Pai et al., 2020). In the atmosphere, OA are either directly emitted as particles from, e.g., volcanic eruption and biomass burning, or formed from gas-to-particle conversion from the oxidation of volatile organic compounds (VOCs). Due to the complex mixture of diverse organic compounds in OA particles, characterization of their chemical composition remains challenging, notably at high time-resolution. Such characterization is crucial for understanding particle formation, growth, aging, as well as their physicochemical properties in the atmosphere. Therefore, improving online detection techniques for OA particles is essential.

Currently, online mass spectrometry (MS), such as chemical ionization mass spectrometer, is one of the key technologies for measuring gaseous oxygenated organic species. It takes advantage of tracking the evolution of compounds during their formation and phase partitioning. To retrieve particle-phase composition with a mass spectrometer, condensed molecules must be converted to gaseous analytes before being ionized. The key point for getting qualitative real-time particle composition information is the design of the inlet in front of the MS. The current technologies to achieve this conversion can be broadly classified into two categories: those requiring pre-collection of particles and those that do not. A thermal desorption chemical ionization mass spectrometer (TD-CIMS) (Voisin et al., 2003; Smith et al., 2004; Li et al., 2021) uses a metal filament in an electrostatic precipitator to collect pre-charged aerosol particles and thereafter thermally evaporate them by pulsing a known current on the

filament. This approach requires both particle deposition and subsequent evaporation. Another instrument involving particle collection is the filter inlet for gases and aerosols (FIGAERO) developed by Lopez-Hilfiker et al. (2014). The FIGAERO inlet collects particles onto a PTFE filter and desorbs them afterward using a heated N₂ flow. This approach enables the control of particle evaporation in a pre-set temperature-ramping program, making it an efficient technique to obtain particle volatility information (Stark et al., 2017; Bannan et al., 2019; Tikkanen et al., 2020; Thornton et al., 2020) while reaching very low detection limits (Lopez-Hilfiker et al., 2014; Thornton et al., 2020). However, since the collection and evaporation of aerosol particles take several tens of minutes, the time resolution of particle measurements remains low, especially in places with low particle loading, resulting in semi-online measurements. Because of particle collection, the FIGAERO is also not interference-free from organic mixtures, especially with high concentrations of aerosol particles (Bannan et al., 2019), which may affect the retrieved volatility of single components. Finally, it has been shown that during evaporation, chemical processes can occur, altering the chemical composition and the information retrieved (Stark et al., 2017; Schobesberger et al., 2018; Buchholz et al., 2020).

In contrast, there are inlets without the need to pre-concentrate particles. An inlet designed for chemical analysis of aerosols online (i.e., CHARON) has been developed without the need to collect particles (Eichler et al., 2015). It consists of a carbon strip denuder for gas-phase compound removal, an aerodynamic lens for particle collimation, and a thermo-desorption unit for particle evaporation. The evaporated compounds can be analyzed with a downstream low-pressure MS. However, due to the strong electric field in the ion drift tube, the protonation-induced ionic fragmentation of oxygenated organic compounds biases the real distribution of particle-phase chemical composition (Müller et al., 2017; Li et al., 2022; Peng et al., 2023). Another technique, using an atmospheric pressure chemical ionization Orbitrap mass spectrometer (APCI), evaporates aerosol particles in a heated ceramic tube where thermal

decomposition compounds can be observed (Vogel et al., 2016; Zuth et al., 2018). Finally, the most recently developed vaporization inlet for aerosols (VIA) coupled to a NO₃-CIMS (VIA-NO₃-CIMS) allows continuous thermal desorption and online detection of particle-phase highly oxidized molecules without pre-concentration (Häkkinen et al., 2023), which mainly consists of a sulfonert-coated stainless steel as TD unit and a following cold dilution flow of N₂. The evaporation tube of VIA is bonded with an insert silica layer into the surface; a dilution flow is used to cool down the sampling flow and minimize the recondensation of the evaporated compounds. The parameters of the dilution unit are critical factors that affect the final sensitivity of the entire system. While this new coupling allows the identification of highly oxygenated organic molecules (HOMs) at atmospheric relevant particle concentrations, the design of the VIA yields subsequent thermal fragmentation when the analytes interact with the heated walls of the TD unit (Zhao et al., 2023; Zhao et al., 2024b). Currently, there is no technique based on thermal evaporation able to prevent thermal fragmentations, and suitable for the on-line measurement of moderate oxygenated (e.g., molecular oxygen atoms <6) organic species.

To prevent thermal fragmentation, the extractive electrospray ionization time-of-flight mass spectrometer (EESI-TOF-MS) (Lopez-Hilfiker et al., 2019) was developed as an online method for particle analysis without the need of thermal desorption. In the EESI-TOF-MS, sampled particles collide with charged electrospray droplets, and the soluble compounds are extracted and ionized through adduct formation. However, the quantification of molecules remains challenging due to the uncertainties in the dependence of instrument sensitivity on molecular identity (Lopez-Hilfiker et al., 2019; Wang et al., 2021). This method exhibits also important background due to the low selectivity of the reagent ions (e.g., Na⁺), making the identification and the quantification of the compounds of interest challenging (Lee et al., 2020; Bell et al., 2023). Additionally, for highly oxygenated and multifunctional compounds,

volatility determination remains particularly uncertain, as isomerism and intermolecular interactions can significantly influence evaporation behavior (Lee et al., 2014; Bannan et al., 2019).

In this work, we designed a newly designed wall-free particle evaporator (WALL-E) to on-line evaporate organic particles while preventing ionization-induced fragmentation and minimizing thermal decomposition effects. WALL-E is coupled to a chemical ionization inlet attached to a CIMS to achieve real-time measurements of aerosol particles at molecular level. An extensive characterization of the WALL-E system is presented here, where its performance is tested as a function of various parameters, e.g., flow rates and evaporation temperatures. The sensitivity of WALL-E coupled to an atmospheric pressure CIMS is determined and the system is used to quantitatively retrieve the composition of particle-phase oxygenated molecules generated from the O₃/OH initiated oxidation of α -pinene in an aerosol flow tube reactor. Finally, by scanning the WALL-E temperature, volatility information can be extracted that can be inferred from the measured thermograms. Polyethylene glycol (PEG) is used to evaluate the volatility measurements and compare our results to the existing techniques.

2 WALL-E design, Experiments, and Instruments

2.1 WALL-E setup and simulations

The WALL-E system is designed to thermally desorb aerosol particles while minimizing the analyte wall-interactions with the system, thus limiting fragmentation (Figure 1). The system enables real-time mass spectrometric analysis of aerosol particles by integrating a series of carefully designed components, including a gas-phase denuder, a thermal desorber (TD) unit equipped with a sheath flow, a ceramic spacer for thermal isolation, and a dilution/cooling unit. Each component is optimized to enhance sample stability and compatibility with different types of CIMS.




Figure 1: (A) Design of the WALL-E interface. (B) Schematic of WALL-E with a gas-phase denuder (GPD) connected to the inlet. The thermal desorption region (TD) is shown in red, where the hot flow (HF) is mixed with the sample flow (SF). A ceramic spacer is indicated by the dashed rectangle. The cooling region (CR), shown in blue, is where the cooling flow (CF) is introduced.

2.1.1 Gas-phase denuder

Like the EESI or the VIA design, the gas-phase denuder (GPD) is the first stage of the WALL-E system (Figure 1B), designed to selectively remove VOCs and other inorganic gaseous species (e.g., nitric acid) while allowing aerosol particles to pass through with minimal losses (transmission efficiency > 90%). This denuder consists of a 10 mm outer diameter, 6 mm inner diameter, stainless steel tube lined with an activated charcoal honeycomb structure, which provides a high surface area for adsorption. Maintaining laminar flow is critical for the denuder's efficiency. With a sample flow rate of 1 SLPM (standard liters per minute), the flow

remains laminar with a Reynolds number of 239, ensuring optimal gas-phase diffusion to the adsorbent walls.

2.1.2 Core vaporization unit

The TD is the core of WALL-E (Figure 1B), where aerosol particles are converted into gas-phase species through flash evaporation. It comprises a 4 cm stainless steel tube heated to up to 390°C. As the sample flow enters the TD region, it merges with a hot nitrogen (N₂, up to 390°C) sheath flow HF (Figure 1B), introduced upstream at a flow rate ranging from 0 to 1 SLPM. All heated regions are insulated to minimize heat exchange with the ambient environment. The hot sheath-flow serves multiple purposes: it maintains a laminar flow by preventing turbulences due to T differences of the sampling flow with the walls of the TD and protects the evaporated analytes from the heated walls. Hence, the TD's design prioritizes minimal wall interactions to reduce potential fragmentation of the products, while ensuring that particle evaporation occurs within a well-defined thermal environment. The stainless-steel tube provides consistent heat transfer along its length, creating a stable thermal zone for the analytes. A limitation of the current design is the use of uncoated stainless steel, which can lead to adsorption or memory effects for more-volatile compounds. Inert coatings will be considered in future iterations to further minimize wall interactions. The combination of the heated tube and the hot sheath-flow allows the aerosol particles to rapidly reach the target temperature, significantly reducing the residence-time required for complete desorption. This rapid temperature-ramp is essential for achieving flash evaporation, where particles are quickly vaporized without prolonged exposure to high temperatures that could lead to thermal decomposition. At a sample flow rate of 1 SLPM, the residence time in the TD is estimated to be between 70 and 90 ms, depending on the added hot dilution flow (0 to 0.5 SLPM).

Numerical simulations using COMSOL Multiphysics® are used to evaluate and refine the flow and temperature dynamics within the system, employing models for turbulent flow,

compressible flow, heat transfer, and chemical transport. These simulations show the system's ability to maintain laminar flow conditions across all units, which is essential to achieve uniform heating and consistent particle transport. Experimental validation is performed in conjunction with these simulations, offering valuable insights into both the operational efficiency and areas for future optimization. The numerical simulations indicate the presence of uniform heating along the tube length, with thermal stabilization achieved within the first centimeter (Figure 2, section A). These simulations highlight a well-defined temperature gradient that optimizes the desorption processes by ensuring a uniform thermal environment throughout the TD. The combination of controlled flow rates, precise heating, and minimized wall-interactions results in efficient aerosol evaporation with minimal fragmentation (Figure 2).




Figure 2: Composite simulation results of the WALL-E system. A: Temperature profile in the TD region, showing efficient heat transfer and controlled mixing as well as the laminar flow profile. B:B': Radial temperature distribution at the TD exit, showing uniform heating across

most of the sample flow. C: Short region of turbulences which allows mixing the heated sample flow with the cooling dilution flow.

2.1.3 Dilution/cooling unit

After the flash evaporation, a dilution/cooling flow (CF) is used by introducing nitrogen (room temperature) via the second dilution/cooling region (CR) to prevent re-condensation of vaporized species (Figure 1B). This dilution unit is also a critical step to reduce turbulences caused by temperature gradients between the TD and the downstream regions as previously observed with the VIA-CIMS system (Zhao et al., 2024b). A ceramic spacer is positioned between the TD and the cooling unit to ensure thermal isolation (Figure 1B). This spacer, made from Alumina Ceramic Al_2O_3 with low thermal conductivity, minimizes heat transfer between the heated TD and the cooled sample stream, maintaining distinct thermal zones, as depicted in Figure 2. Simulations confirm the effectiveness of the ceramic spacer in preserving thermal gradients, greatly reducing unwanted heat transfer into the cooling unit. To further optimize the cooling unit, two fans are mounted right after the second dilution region to ensure that the sample reaches room temperature while preventing turbulence within the CI inlet. This separation is critical to ensure the stability of the sample as it progresses toward the CI inlet operated at atmospheric pressure. As indicated in Figure 2 section C, the cooling seems to be achieved within the first centimeter while minimizing turbulence along the tube .

2.2 Experiments

2.2.1 Characterization of the optimal WALL-E setting parameters

To determine the optimum WALL-E parameters (sampling flow rate, hot sheath-flow/temperature, and TD temperature), particles are generated using an atomizer (Model 3076, TSI, Minnesota, USA) with an aqueous solution of citric acid (Sigma-Aldrich, $\geq 99.5\%$), d-(+)-glucose (Alfa Aesar, $\geq 99\%$), malonic acid (Sigma-Aldrich, 99%), phthalic acid (Sigma-Aldrich,

$\geq 99.5\%$), and ammonium sulfate (Sigma-Aldrich, $\geq 99\%$). The concentration of each chemical in the aqueous solution is 1 ppm. Aerosol particles are dried using a silica gel dryer, after which the sampling line is divided to provide an aerosol flow to a scanning mobility particle sizer (SMPS, TSI Incorporated, USA) and to the WALL-E CIMS. The setup is illustrated in Figure S1 and all parameters tested are summarized in Table S1.

2.2.2 Sensitivity determination

The sensitivity of individual compounds is directly proportional to their clustering strength with the reagent ions (Iyer et al., 2016; Bi et al., 2021), which can be probed by performing a declustering scanning procedure (Lopez-Hilfiker et al., 2016). Using the setup shown in Figure S1, single component aerosol particles are generated using an atomizer containing single component aqueous solutions of the following compounds: 1,5-dihydroxynaphthalene (Sigma-Aldrich, $\geq 97\%$), 3,4,5-trihydroxybenzaldehyde (Sigma-Aldrich, $\geq 98\%$), 4-hydroxyphenylacetic acid (Sigma-Aldrich, $\geq 98\%$), ammonium sulfate (Sigma-Aldrich, $\geq 99\%$), citric acid (Sigma-Aldrich, $\geq 99.5\%$), d-(+)-glucose (Alfa Aesar, $\geq 99\%$), d-mannitol (Sigma-Aldrich, $\geq 98\%$), phthalic acid (Sigma-Aldrich, $\geq 99.5\%$), phthalic acid d4 (Sigma-Aldrich, $\geq 98\%$), and shikimic acid (Alfa Aesar, $\geq 98\%$). The sensitivity, reported in Table S2, for each type of particle is determined by the linear regression of the mass concentration measured by the SMPS and the normalized signal intensity of the analyte clustered with Br^- (i.e., $[\text{M-Br}^-]$) detected by the WALL-E-CIMS (Figure S2). To further assess the correlation between the sensitivity and the binding energy of the detected ion adducts, voltage scanning (i.e., increasing the voltage difference between two ion optics) is performed to determine the half-signal maximum intensity (Lopez-Hilfiker et al., 2016; Riva et al., 2020) for the generated single-component aerosol particles.

2.2.3 SOA particle generation

To further examine the performance of the WALL-E system, SOA are generated in an 18-liter Pyrex glass aerosol flow tube reactor (12 cm i.d. \times 158 cm length) from the O₃/OH initiated oxidation of α -pinene in the presence of SO₂ at room temperature and atmospheric pressure (Stein and Scott, 1994). Ozone (O₃) is stably generated by passing a flow of 0.6 SLPM of synthetic air after exposure to a UV lamp (Ozone Generator Model 610, Jelight Company, Inc, Irvine, USA). SO₂ is injected from a commercial cylinder (500 ppm, AIR PRODUCTS Inc.), to promote the particle formation and to generate more SOA mass. α -Pinene was introduced from a pressurized cylinder (40 ppm in nitrogen). O₃, SO₂, and α -pinene are continuously injected into the aerosol flow tube reactor to generate a total aerosol mass ranging from 1.0 to 15.6 $\mu\text{g}\cdot\text{m}^{-3}$. The concentrations of reactants are summarized in Table S3. A mixture of nitrogen and oxygen (total flow 4 SLPM) is used as a carrier gas, providing a reaction time of \sim 4.5 minutes. Before injecting α -pinene, background measurements are obtained.

2.2.4 Thermograms and T_{max} determination

The volatility distribution of aerosol particles is investigated using thermograms obtained with the WALL-E system. Polyethylene glycol (PEG-400) aerosols are produced by atomizing aqueous solutions as performed in previous studies (Ylisirniö et al., 2021; Zhao et al., 2024b) and sampled by the WALL-E system. The temperature of the WALL-E TD is gradually increased from 30°C to 390°C in 30°C increments every 10 minutes. The PEG standards, ranging from PEG-6 to PEG-17, have been chosen to represent a wide range of molecular weights and volatilities (Krieger et al., 2018).

2.3 Instrumentation

Within these experiments, WALL-E is associated with an atmospheric pressure CI inlet (Riva et al., 2019a; Riva et al., 2020) coupled to an Orbitrap (Q-Exactive, Thermo Fisher Scientific)

utilizing bromide ions (Br^-) as the reagent ion. Br^- is generated from dibromomethane (Sigma-Aldrich, 99%) continuously flushed by 2 standard cubic centimeters per minute (sccm) of pure N_2 , and subsequently ionized with a soft X-ray photoionizer (Hamamatsu, L9491). The sheath and the total flows are 24 SLPM and 33.5 SLPM, respectively. The Orbitrap is operated with an automatic gain control (AGC target) of 1×10^6 charges, an S-lens radio frequency level of 50, a maximum injection time of 1000 ms, 10 microscans, and a capillary temperature set to 150 °C. The mass resolution is 140,000 (at m/Q 200). Orbitool 2.2.4 (Cai et al., 2021) is used for analyzing the data. The data are pre-averaged to 1 minute. Signals are background subtracted and normalized by the signal intensity of Br^- (m/Q 79). To obtain an accurate concentration of compounds present in low abundance, a linearity correction (Riva et al., 2020) is applied to all measured signals (Figure S3). All product ions are normalized to the Br^- signals.

The mass concentration of particles is retrieved using a scanning mobility particle sizer (SMPS) utilizing a differential mobility analyzer (DMA; 3081, TSI Inc.) connected to a CPC (3772, TSI Inc.), by applying an assumed particle density (1.45 g/cm^3) (Kim et al., 2010; Shilling et al., 2009) for aerosol particles generated from the O_3/OH initiated oxidation of α -pinene in the presence of SO_2 . An impactor (0.071 m) is used, and the sampling flow rate is 1 SLPM with a sheath flow of 10 SLPM.


3 Results and Discussion

3.1 Optimal setting parameters of WALL-E

As described in section 2.1, WALL-E involves multiple operational parameters, including the sample flow rate (SF), the hot sheath flow rate (HF) and its temperature upstream of the TD, the dilution flow in the cooling region, and the TD temperature region. The cooling flow is held constant at 10 SLPM to minimize turbulence in the CI inlet operated at atmospheric pressure. To determine the optimized parameters, various combinations of SF and HF flowrates as well

as HF and TD temperatures are tested on aerosol particles generated from an atomized water solution containing a mix of standards. All data are corrected for their respective dilution factors. Figure 3 shows the summed signal of all evaporated products (i.e., phthalic acid, citric acid, malonic acid, glucose, sulfuric acid) under these different conditions, normalized to the maximum value.

Figure 3: (A) Relative signal at a fixed hot sheath-flow (HF) temperature of 300 °C versus TD temperature for various combinations of sample flow rate (SF) and HF. (B) Relative signal at a fixed TD temperature of 300 °C versus HF temperature. (C) Relative signal versus SF at



different HF settings. (D) Relative signal versus HF at different SF values. Both (C) and (D) are at a fixed TD and HF temperature of 300 °C. SF is the sample flow rate delivered to the TD, while HF is the heated dilution flow added upstream of the ambient dilution stage. The second, ambient dilution flow and the downstream cooling flow remain fixed to minimize turbulence in the CIMS. Error bars represent standard deviations based on repeated measurements at the same conditions, and all values are corrected for their respective dilution factors.

As depicted in Figure 3A, it is evident that a SF of 1.0 L min⁻¹ consistently results in higher signals than 0.5 L min⁻¹ across all TD temperatures. A SF of 1.5 L min⁻¹ also performs reasonably well, whereas the 0.5 L min⁻¹ configurations exhibit significantly lower signals. A similar trend is observed in Figure 3B, where the relative signal is plotted against the HF temperature. The results indicate that efficient evaporation occurs at a TD temperature of around 300-350°C, while variations in the HF temperature have a less important effect. However, the presence of a HF (either 0.25 or 0.5 L min⁻¹) at an SF of 1 or 1.5 L min⁻¹ enhances evaporation efficiency and flow stability. This influence of flow conditions is more distinctly observed in Figures 3C and 3D and further illustrates the impact of SF and HF on the relative signal. Figure 3C presents the relative signal as a function of SF for different HF settings, demonstrating that a SF of 1 L min⁻¹ and a HF of 0.25 L min⁻¹ yield higher signal intensities compared to other conditions. Similarly, Figure 3D shows the relative signal plotted against HF at different SF values, reinforcing this observation.

A key consideration when selecting optimal SF and HF conditions is to keep fragmentations minimal. Excessive heating in the TD can enhance fragmentation, potentially resulting in thermal decomposition of the analytes. To assess this, we evaluate the thermal decomposition of citric acid, which is a known analyte to decompose within TD inlets (Yang et al., 2021). As shown in Figure S4, the thermal fragmentation of citric acid using the conditions (i.e., SF of 1 L min⁻¹ and a HF of 0.25 L min⁻¹) described above remains negligible, with up to 2% of the total signal attributed to fragment ions at the hottest temperature (i.e., 390°C). This confirms the minimal fragmentation while keeping an effective evaporation.

3.2 Characterization of SOA derived from α -pinene ozonolysis

3.2.1 Raw mass spectra

To underline the performance of WALL-E in measuring and quantifying a complex mixture of OA, the characterization of SOA generated from the oxidation of α -pinene in the presence of SO_2 is used as an example. Figures 4 and S5 present the mass spectra of α -pinene-derived SOA with particle mass concentrations ranging from 1.0 to 15.6 $\mu\text{g m}^{-3}$. These mass concentrations cover atmospherically relevant ranges observed in remote and rural places (Jimenez et al., 2009), which is achieved when using the optimal setting parameters determined in section 3.1. A total of 146 organic monomers ($\text{C}_{1-10}\text{H}_y\text{O}_z\text{S}_0$) and 206 dimers ($\text{C}_{11-20}\text{H}_y\text{O}_z\text{S}_0$) are identified, respectively contributing 67-74 % and 29-23 % to the total signal intensity. The dominant ions identified in the monomer region are assigned to $\text{C}_{8-10}\text{H}_{12-18}\text{O}_{3-7}$ compounds with high-resolution peak fitting in the range of m/Q at 249-343 Th, and $\text{C}_{17-20}\text{H}_{26-34}\text{O}_{5-11}$ compounds in the dimeric range between 405-525 Th. Among the monomers, $\text{C}_{10}\text{H}_{14,16}\text{O}_{3-7}$ are the most abundant ones, followed by $\text{C}_9\text{H}_{14}\text{O}_{3-7}$. The product distribution measured is consistent with previous studies (Zhang et al., 2015; Kahnt et al., 2018; Zhao et al., 2023; Zhao et al., 2024a). It should be mentioned that $\text{C}_8\text{H}_{12}\text{O}_4$ is the most abundant individual compound, which could be norpinic acid or terpenylic acid identified by previous studies (Zhang et al., 2015; Du et al., 2022; Witkowski et al., 2023). It has an unexpected high background of ~ 0.04 normalized counts per second (ncps) when there is no VOC injected. Desorption from the walls of the aerosol flow tube reactor when aerosol particles are produced can explain the presence of $\text{C}_8\text{H}_{12}\text{O}_4$ as reported in previous laboratory studies (Riva et al., 2019b; Wong et al., 2022).

Without sensitivity correction, compounds are significantly detected at a total signal of 0.04 ncps when the particle mass is $1.0 \pm 0.1 \mu\text{g m}^{-3}$. As particle mass concentration increases, the normalized signals of particle-phase products also rise, reaching 0.8 ncps at $15.7 \pm 1.2 \mu\text{g m}^{-3}$. The detected products exhibit a strong linear response to particle mass concentration, as

shown in Figure S6 illustrating the suitability of the WALL-E system at measuring aerosol particles under atmospheric relevant conditions.




Figure 4. Mass spectra of particle-phase organic compounds formed from the oxidation of α -pinene under varying mass particle concentrations (1.0 - $15.6 \mu\text{g}\cdot\text{m}^{-3}$). Compounds are evaporated and detected by the WALL-E CIMS with HF and TD temperatures of 300°C , SF of 1 SLPM , HF of 0.25 SLPM , and a cold dilution flow of 10 SLPM . Signals are normalized to Br^- signal intensity. Red and blue refer to monomers ($\text{C}_{1-10}\text{H}_y\text{O}_x$) and dimers ($\text{C}_{11-20}\text{H}_y\text{O}_x$), respectively. The normalized signals of dimeric compounds are multiplied by a factor of 10.

3.2.2 Sensitivity determination on SOA molecules

The sensitivity of standard compounds, which is determined by the linear regression between the normalized signals detected by WALL-E and the particle mass concentrations (Figure S2), is summarized in Table S2. The r^2 -values of the fitting for most standard compounds are good (0.93-0.99), while 1,5-dihydroxy naphthalene has the lowest r^2 -value of 0.85. By utilizing the in-source collision ion dissociation feature (Riva et al., 2019a), which corresponds to an increase in the direct current (DC) offset voltages between two ion optics within the flatapole, the binding energy of the $[M-Br^-]$ adducts can be probed (Figure S7). Their voltage difference at half signal maximum intensity (dV_{50}) broadly ranges from 5.1 to 20.2 Volts, indicating differences in binding energies and varying clustering strengths to Br^- . Consequently, the correlation between the sensitivity and the dV_{50} values obtained from the standard compounds is fitted by a non-linear sigmoidal function (Figure S8), which is consistent with prior studies using the same approach to quantify gaseous species (Lopez-Hilfiker et al., 2016; Iyer et al., 2016; Zaytsev et al., 2019; Xu et al., 2022). This calibration curve provides an estimation of the system sensitivity based on experimentally obtained dV_{50} cluster values from declustering scans with increasing energy (Lopez-Hilfiker et al., 2019). Using this method, raw MS signal intensities can be converted into quantified amounts, reducing the need for compound-specific calibration when authentic standards are not available. This method enables semi-quantification across a wide variety of molecules.

The correlation between sensitivity (Figure S2) and dV_{50} (Figure S7) is obtained based on the standard compounds, as depicted in Figure S8. By applying this sigmoidal function and using the dV_{50} values determined for individual α -pinene-derived SOA products, the concentration of every oxidation compound can be estimated. As an example, Figure S9 presents the declustering profiles of C_{10} monomers and C_{19-20} dimers from SOA. Consistent with previous studies (Riva et al., 2019a), more oxidized compounds exhibit stronger binding

energies, resulting in higher dV_{50} values. As shown in Figure S10, the corresponding sensitivity of α -pinene-derived SOA compounds generally increases with molecular mass and reaches a plateau corresponding to the maximum sensitivity (i.e., collision limit), with an upper limit of sensitivity of $0.08 \text{ ncps}\cdot\text{per}\cdot\mu\text{g}\cdot\text{m}^{-3}$ providing a limit of quantification $\sim 5\text{-}10 \text{ pg}\cdot\text{m}^{-3}$ corresponding to uncorrected signal intensities of $5 \times 10^{-6} - 1 \times 10^{-5}$ (Riva et al., 2020).

3.2.3 Sensitivity-corrected chemical composition of SOA particles

By assessing the sensitivity of individual α -pinene-derived SOA compounds, the mass concentrations of all identified particle-phase oxidation products can be estimated (Figure S11). Using the correlation between sensitivity and dV_{50} based on 10 standard compounds, the total particle mass concentration is estimated to be $27.1 \mu\text{g}\cdot\text{m}^{-3}$, 74% higher than the SMPS measurements ($15.6 \mu\text{g}\cdot\text{m}^{-3}$), assuming spherical particles with an aerosol density of 1.45 g/m^3 (Kim et al., 2010; Shilling et al., 2009). Among all calibrated compounds, two outliers exist (i.e., shikimic acid and glucose), which might be due to partial evaporation leading to an underestimation of the sensitivity. When excluding these two species, the estimate of total particle mass concentration is closer (with 36% overestimation) to the total particle mass concentration measured by the SMPS as depicted in Figure 5. As discussed in previous studies, selected standard compounds might induce uncertainty in the sensitivity estimations (Zaytsev et al., 2019; Bi et al., 2021; Song et al., 2024). Notably, the presence of different isomers can yield substantial uncertainties, especially when their sensitivity may vary by an order of magnitude (e.g., Lee et al., 2014). It should also be mentioned that the total mass concentration determined by the SMPS is prone to uncertainties (Wilson et al., 2015; Bell et al., 2023). Using the dV_{50} method presented here and considering a total SOA mass of $15.6 \mu\text{g m}^{-3}$, the mass concentration of H_2SO_4 is estimated at $8.7 \mu\text{g m}^{-3}$ (using 10-compound fitting) and $6.4 \mu\text{g m}^{-3}$ (using 8-compound fitting). A good agreement between the direct calibration (Table S2) and the declustering scan method is retrieved (7% overestimation and 25% underestimation for the

10 and 8-compounds fit, respectively underlining the benefit of using this approach to obtain the concentrations of organic and inorganic present in the particles. Overall, a linear relationship exists with the total SOA mass concentrations in the range of 1.0-15.6 $\mu\text{g}\cdot\text{m}^{-3}$ (Figure 5) and exhibits an overall good agreement between WALL-E and SMPS measurements (i.e., 20-30%). Hence, this method can provide a deeper understanding of aerosol composition and evolution. For example, due to their general higher sensitivities, dimers contribute only 14-18% to the total mass (Figures 5, S10), much less than their fractions based on signal intensity (23-29%). This suggests that the dimer contributions to particles may be overestimated when solely based on signal intensity.




Figure 5. (A) The estimated mass concentration measured by Wall-E as a function of SOA mass concentration measured by the SMPS. (B) Fractions of monomers, dimers and sulphur-containing compounds in the SOA particles formed from O_3/OH initiated oxidation of α -pinene in the presence of SO_2 . The calculation of fractions is weighted by normalized signals (dashed bars) and mass concentrations (solid bars), respectively.

3.3 Assessing particle molecular volatility

3.3.1 Thermograms and T_{max} determination

The determination of volatility represents one of the greatest analytical challenges when characterizing aerosol particles, as it depends on multiple factors, including molecular composition, intermolecular interactions, and experimental conditions (Compernelle et al., 2011). Various experimental and theoretical techniques have been developed over the last decades to retrieve volatility information, each with advantages and limitations.

While FIGAERO and VIA have been widely used, their design constraints introduce inherent limitations. FIGAERO can introduce artifacts such as recondensation, analyte interactions, and fragmentation. The prolonged residence time on the FIGAERO may also lead to early desorption of volatile species or chemical reactions between co-deposited compounds, impacting the accuracy of volatility estimates (Stark et al., 2017; Schobesberger et al., 2018; Buchholz et al., 2020). VIA thermograms, on the other hand, show evidence of fragmentation and thermal decomposition at high temperatures (Zhao et al., 2024b). WALL-E introduces a new approach, optimizing the balance between thermal residence time and evaporation efficiency, allowing for precise volatility determination with reduced wall interactions.

3.3.2 T_{50} as a Robust Volatility Metric

The T_{\max} values, which correspond to the peak desorption temperature, represent the temperature at which the maximum desorption rate occurs. Traditionally, T_{\max} has been used to estimate volatility from thermograms. This approach is well-suited for FIGAERO and VIA, where the thermograms typically exhibit a near-Gaussian profile, with a clear peak followed by a steep signal decline due to decomposition and fragmentation. However, in the case of WALL-E, where thermograms more closely follow a sigmoidal desorption trend, T_{\max} determination becomes less straightforward. At high temperatures, the signal does not drop sharply but instead asymptotically reaches a plateau, making T_{\max} highly sensitive to noise and minor variations in the upper temperature range.

To overcome this limitation, we propose using T_{50} , the temperature at which the signal reaches 50% of its maximum as a more robust volatility metric. Since T_{50} is located in the steepest region of the sigmoid fit, it is significantly less affected by noise, baseline shifts, or small variations in signal intensity. Unlike T_{\max} , which depends on the choice of an arbitrary threshold (e.g., 99.5% or 99.9% of the signal maximum), T_{50} is a main feature of the sigmoid function, making it a more consistent and reproducible parameter for comparing volatility trends in the case of WALL-E. To accurately determine T_{50} , the thermograms are smoothed and fitted using a sigmoid function. This approach minimizes the influence of temperature ramping increments and instrumental noise. Figure 6 displays thermograms for PEG standards, illustrating that less volatile compounds require higher temperatures for complete evaporation, whereas more volatile compounds evaporate at lower temperatures.




Figure 6: Thermograms of PEG standards obtained using WALL-E. The desorption profiles illustrate the relative signal intensity as a function of temperature.

In contrast to earlier systems, including the FIGAERO and VIA, which often produce Gaussian-shaped thermograms with clearly defined peaks, WALL-E thermograms follow a sigmoidal profile. This makes the determination of T_{\max} more sensitive to the selected threshold and prone to variability, especially at upper temperature ranges.

The use of a sigmoid fit for the WALL-E data is justified by the observed thermograms, which closely resemble this behaviour, with only a moderate signal decrease at the highest tested temperatures. For comparison purposes with other particle evaporators, T_{\max} is used with a signal representation of 98%. Figure S12 shows the effect of varying this signal percentage threshold. By varying between 99.5% and 98%, a difference of about 10% in the T_{\max} value is observed, which makes it challenging to select an appropriate threshold. In addition, at higher temperatures, the fit for less volatile PEGs can result in more errors due to fewer data points, which mainly affect PEG-17. Taking a threshold around where the signal is half always results in more accurate results even with 10% variations in the threshold. The signals of more volatile compounds, such as PEG-6 and PEG-7, with T_{\max} values of 203.6°C and 246.3°C, decrease by only 11.2% and 18.9% respectively, even after a temperature increase of over 150°C, which represents a significant improvement compared to other online TD techniques. At this point, PEG-16 and PEG-17 reach their T_{\max} (Table S4). This demonstrates WALL-E's ability to maintain signal integrity across a wide volatility range, ensuring reliable thermal desorption without excessive signal loss.

3.3.3 Temperature Correction and Comparison with VIA and FIGAERO

Although T_{50} is introduced as a more robust metric in this study, temperature correction and inter-comparison with other techniques have traditionally been performed using T_{\max} values. Therefore, we present the correction and comparison based on T_{\max} to ensure consistency with previous studies before applying the T_{50} approach.

The temperatures applied during thermal desorption in WALL-E do not directly reflect the actual gas-phase temperatures experienced by the desorbed compounds. Due to thermal lag and heat transfer dynamics, a correction factor was derived based on the COMSOL CFD simulations, as well as actual measurements using a thermocouple at the end of the TD. The internal gas temperature at the core of the TD was found to follow a linear relationship with the


set temperature (Figure S13). This correction ensures that the reported T_{\max} values accurately represent the true volatility behaviour of the analytes. The corrected and uncorrected T_{\max} and T_{50} values for each PEG compound are plotted against their molecular masses in Figure S14. The corrected values align well with those reported for VIA and FIGAERO, maintaining a consistent trend across different molecular weights (Figure S14). In contrast, the uncorrected T_{\max} values are systematically higher, displaying a different slope compared to the other two systems. This highlights the necessity of applying the correction factor to ensure accurate volatility determination.

T_{\max} values from WALL-E and VIA are higher compared to FIGAERO, which can be attributed to the short residence time in the TD region (Zhao et al., 2024b). Additionally, WALL-E exhibits lower T_{\max} values than VIA, likely due to the introduction of heated dilution flow into the sample stream. This promotes flash evaporation, causing compounds to desorb before reaching full thermal equilibrium, ultimately shifting T_{\max} to lower values while maintaining the expected volatility trend. As illustrated in Figure S15, the residence time in the TD region further influences T_{50} or T_{\max} values. Longer residence times allow for gradual heating and equilibration, leading to lower temperature values. In contrast, shorter residence times accelerate desorption, resulting in higher temperature values due to insufficient thermal equilibration.

3.3.4 Volatility Estimation from T_{50}

As previously explained, the use of T_{50} would provide a more reliable estimation of the volatility, in the case of WALL-E, which is typically inferred using the relationship between T_{50} and the saturation concentration (C^*) as discussed in prior studies (Ylisirniö et al., 2021). The wide range of T_{50} values observed for PEG standards underscores the broad applicability of WALL-E for volatility characterization across diverse classes of compounds. To derive volatility

estimates, we apply the parameterization method proposed previously (Krieger et al., 2018; Ylisirniö et al., 2021), which utilizes measured vapor pressures for PEG-5 to PEG-8 and




extrapolates the trend for higher molecular weights. Additionally, an alternative approach using the parameterization (Li et al., 2016) is considered. While both methods produce similar trends for lower-mass PEGs, they diverge significantly at higher masses, reflecting the inherent uncertainties in extrapolating volatility predictions. Given these discrepancies, we define a volatility range (Figure 7) that encompasses both parameterizations, providing a more robust estimation framework until additional direct vapor pressure measurements become available.

Figure 7: Volatility estimation using measured T_{50} values, ranging from SVOC to ULVOC.

To assess the applicability of T_{50} -based volatility determination, we apply this approach to α -pinene derived SOA compounds, specifically $C_{10}H_{16}O_{6-9}$. The T_{50} values of these SOA species are determined using multiple ramp-up speeds, ensuring reproducibility across different heating rates. Additionally, a faster temperature ramp-down from 390°C to 30°C is used as an independent validation method. As shown in Figure 8, the T_{50} estimated using WALL-E for the

$C_{10}H_{16}O_{6-9}$ from the heating and cooling phases (i.e., temperature ramping up and down, respectively) exhibit a very good agreement confirming the robustness of the system.

Determined T_{50} values are further compared with theoretical volatility estimation models, specifically COSMO-RS and SIMPOL, as presented previously (Kurtén et al., 2016; Peräkylä et al., 2020). Our measured T_{50} values for the SOA compounds fall well within the SIMPOL-predicted region, confirming that the volatility estimates obtained using WALL-E are consistent with theoretical predictions. As mentioned earlier, without a pre-separation method,



the presence of isomers can alter the quantification of the compounds of interest, which is also the case for volatility estimation. As previously discussed by Kurtén et al. (2016) and Peräkylä et al. (2020), isomeric structures can exhibit significantly different volatilities, reinforcing the need to account for molecular configurations beyond elemental composition when interpreting volatility trends.

Figure 8: Comparison of estimated T_{50} values for α -pinene derived SOA compounds, specifically $C_{10}H_{16}O_{6-9}$ with volatility predictions from COSMO-RS and SIMPOL models. The measured T_{50} values (orange) are shown with their min/max range, while the shaded regions

represent model predictions taking into consideration the isomerisation. We applied the sigmoid fitting method to both the heating, cooling, and fast cooling ramps, extracted the T_{50} values from each, and averaged them. The orange line represents this averaged T_{50} value.

4 Conclusions

Within this work we present a new analytical method (i.e., WALL-E) to retrieve the chemical composition of atmospheric particles in real time. Coupled with a CIMS using Br^- ion chemistry as the reagent ion, WALL-E is comprehensively characterized, achieving efficient particle evaporation with maximum evaporation efficiency while exhibiting minimal thermal decomposition across a range of operational settings. The characterization of SOA produced from the O_3/OH initiated oxidation of α -pinene in the presence of SO_2 further validates WALL-E's performance in resolving and quantifying complex organic aerosol mixtures. It successfully identifies a broad range of monomeric (C_{1-10}) and dimeric (C_{11-20}) compounds. By utilizing the in-source collision ion dissociation feature, the sensitivity of the analytical method is realized using a wide variety of authentic standards used to determine the correlation between binding energy and sensitivity. By using this function and performing declustering procedures, individual α -pinene derived SOA compounds are quantified at concentration as low as $10 \text{ pg} \cdot \text{m}^{-3}$ (for a total SOA mass of $1 \text{ } \mu\text{g} \cdot \text{m}^{-3}$). The total estimated SOA mass concentrations is in good agreement with particle concentration measurements obtained by an SMPS, demonstrating the benefit of this approach. Notably, the mass contribution of dimeric compounds is determined, which reveals that they account for only 14-18% of total particle mass, which is notably lower than their fractions (23-29%) based on signal intensity. Finally, the volatility assessment using thermogram analysis demonstrates WALL-E's capability to retrieve T_{50} values with high precision, aligning well with predicted SIMPOL volatility. Future studies should focus on systematically characterizing SOA volatilities across a broader range of precursor compounds

and oxidation conditions, leveraging complementary mass spectrometry and computational modelling techniques to refine volatility estimation approaches.

Overall, WALL-E represents a useful and promising tool for atmospheric research, bridging important gaps in real-time aerosol characterization, quantification of chemical composition of complex particulate organic mixtures, and volatility assessment. It improves the time resolution and minimizes measurement artifacts notably due to thermal fragmentation, providing a new technique for investigating the real-time changes in the formation and growth of atmospheric particles for laboratory and field observations.

Data availability

Data of all figures and tables are available on request from the corresponding author.

Author contributions

MR, SP, IZ, MD, and FB designed and built the WALL-E. LG, IZ, ES, CC, and FSD conducted the experiments. LG, IZ, ES, and FSD analyzed the data. IZ did the simulation. LG and IZ prepared the paper with contributions from all co-authors.

Competing interests

The authors declare no competing financial interest.

Acknowledgement

The authors thank Dr. Siegfried Schobesberger for discussing and sharing the FIGAERO measurements and Dr. Georgios Gkatzelis for discussing the declustering scan.

Financial support

596 This work is funded by the European Research Council Grant (ERC-StG MAARvEL; 423 No.
597 852161). IZ acknowledges funding from the CLOUD-DOC project (Grant Agreement No.
598 101073026) under the HORIZON-MSCA-2021-DN-01 program.

599

600 References

- 601 Bannan, T. J., Le Breton, M., Priestley, M., Worrall, S. D., Bacak, A., Marsden, N. A., Mehra, A., Hammes,
602 J., Hallquist, M., Alfarra, M. R., Krieger, U. K., Reid, J. P., Jayne, J., Robinson, W., McFiggans, G., Coe, H.,
603 Percival, C. J., and Topping, D.: A method for extracting calibrated volatility information from the
604 FIGAERO-HR-ToF-CIMS and its experimental application, *Atmos. Meas. Tech.*, 12, 1429-1439,
605 10.5194/amt-12-1429-2019, 2019.
- 606 Bell, D. M., Zhang, J., Top, J., Bogler, S., Surdu, M., Slowik, J. G., Prevot, A. S. H., and El Haddad, I.:
607 Sensitivity Constraints of Extractive Electrospray for a Model System and Secondary Organic Aerosol,
608 *Analytical Chemistry*, 95, 13788-13795, 10.1021/acs.analchem.3c00441, 2023.
- 609 Bi, C., Krechmer, J. E., Canagaratna, M. R., and Isaacman-VanWertz, G.: Correcting bias in log-linear
610 instrument calibrations in the context of chemical ionization mass spectrometry, *Atmos. Meas. Tech.*,
611 14, 6551-6560, 10.5194/amt-14-6551-2021, 2021.
- 612 Buchholz, A., Ylisirniö, A., Huang, W., Mohr, C., Canagaratna, M., Worsnop, D. R., Schobesberger, S.,
613 and Virtanen, A.: Deconvolution of FIGAERO-CIMS thermal desorption profiles using positive matrix
614 factorisation to identify chemical and physical processes during particle evaporation, *Atmos. Chem.*
615 *Phys.*, 20, 7693-7716, 10.5194/acp-20-7693-2020, 2020.
- 616 Cai, R., Li, Y., Clément, Y., Li, D., Dubois, C., Fabre, M., Besson, L., Perrier, S., George, C., Ehn, M., Huang,
617 C., Yi, P., Ma, Y., and Riva, M.: Orbitool: a software tool for analyzing online Orbitrap mass spectrometry
618 data, *Atmos. Meas. Tech.*, 14, 2377-2387, 10.5194/amt-14-2377-2021, 2021.
- 619 Charnawskas, J. C., Alpert, P. A., Lambe, A. T., Berkemeier, T., O'Brien, R. E., Massoli, P., Onasch, T. B.,
620 Shiraiwa, M., Moffet, R. C., Gilles, M. K., Davidovits, P., Worsnop, D. R., and Knopf, D. A.: Condensed-
621 phase biogenic–anthropogenic interactions with implications for cold cloud formation, *Faraday*
622 *Discussions*, 200, 165-194, 10.1039/C7FD00010C, 2017.
- 623 Compennolle, S., Ceulemans, K., and Müller, J. F.: EVAPORATION: a new vapour pressure estimation
624 method for organic molecules including non-additivity and intramolecular interactions, *Atmos. Chem.*
625 *Phys.*, 11, 9431-9450, 10.5194/acp-11-9431-2011, 2011.
- 626 Du, M., Voliotis, A., Shao, Y., Wang, Y., Bannan, T. J., Pereira, K. L., Hamilton, J. F., Percival, C. J., Alfarra,
627 M. R., and McFiggans, G.: Combined application of online FIGAERO-CIMS and offline LC-Orbitrap mass
628 spectrometry (MS) to characterize the chemical composition of secondary organic aerosol (SOA) in
629 smog chamber studies, *Atmos. Meas. Tech.*, 15, 4385-4406, 10.5194/amt-15-4385-2022, 2022.
- 630 Eichler, P., Müller, M., D'Anna, B., and Wisthaler, A.: A novel inlet system for online chemical analysis
631 of semi-volatile submicron particulate matter, *Atmos. Meas. Tech.*, 8, 1353-1360, 10.5194/amt-8-
632 1353-2015, 2015.
- 633 Fehsenfeld, F., Calvert, J., Fall, R., Goldan, P., Guenther, A. B., Hewitt, C. N., Lamb, B., Liu, S., Trainer,
634 M., Westberg, H., and Zimmerman, P.: Emissions of volatile organic compounds from vegetation and
635 the implications for atmospheric chemistry, *Global Biogeochemical Cycles*, 6, 389-430,
636 <https://doi.org/10.1029/92GB02125>, 1992.
- 637 Häkkinen, E., Zhao, J., Graeffe, F., Fauré, N., Krechmer, J. E., Worsnop, D., Timonen, H., Ehn, M., and
638 Kangasluoma, J.: Online measurement of highly oxygenated compounds from organic aerosol, *Atmos.*
639 *Meas. Tech.*, 16, 1705-1721, 10.5194/amt-16-1705-2023, 2023.
- 640 Iyer, S., Lopez-Hilfiker, F., Lee, B. H., Thornton, J. A., and Kurtén, T.: Modeling the Detection of Organic
641 and Inorganic Compounds Using Iodide-Based Chemical Ionization, *The Journal of Physical Chemistry*
642 *A*, 120, 576-587, 10.1021/acs.jpca.5b09837, 2016.
- 643 Jimenez, J. L., Canagaratna, M. R., Donahue, N. M., Prevot, A. S. H., Zhang, Q., Kroll, J. H., DeCarlo, P.
644 F., Allan, J. D., Coe, H., Ng, N. L., Aiken, A. C., Docherty, K. S., Ulbrich, I. M., Grieshop, A. P., Robinson,
645 A. L., Duplissy, J., Smith, J. D., Wilson, K. R., Lanz, V. A., Hueglin, C., Sun, Y. L., Tian, J., Laaksonen, A.,
646 Raatikainen, T., Rautiainen, J., Vaattovaara, P., Ehn, M., Kulmala, M., Tomlinson, J. M., Collins, D. R.,
647 Cubison, M. J., E., Dunlea, J., Huffman, J. A., Onasch, T. B., Alfarra, M. R., Williams, P. I., Bower, K.,
648 Kondo, Y., Schneider, J., Drewnick, F., Borrmann, S., Weimer, S., Demerjian, K., Salcedo, D., Cottrell, L.,
649 Griffin, R., Takami, A., Miyoshi, T., Hatakeyama, S., Shimono, A., Sun, J. Y., Zhang, Y. M., Dzepina, K.,
650 Kimmel, J. R., Sueper, D., Jayne, J. T., Herndon, S. C., Trimborn, A. M., Williams, L. R., Wood, E. C.,

Middlebrook, A. M., Kolb, C. E., Baltensperger, U., and Worsnop, D. R.: Evolution of Organic Aerosols in the Atmosphere, *Science*, 326, 1525-1529, doi:10.1126/science.1180353, 2009.

Kahnt, A., Vermeylen, R., Iinuma, Y., Safi Shalamzari, M., Maenhaut, W., and Claeys, M.: High-molecular-weight esters in α -pinene ozonolysis secondary organic aerosol: structural characterization and mechanistic proposal for their formation from highly oxygenated molecules, *Atmos. Chem. Phys.*, 18, 8453-8467, 10.5194/acp-18-8453-2018, 2018.

Kanakidou, M., Seinfeld, J. H., Pandis, S. N., Barnes, I., Dentener, F. J., Facchini, M. C., Van Dingenen, R., Ervens, B., Nenes, A., Nielsen, C. J., Swietlicki, E., Putaud, J. P., Balkanski, Y., Fuzzi, S., Horth, J., Moortgat, G. K., Winterhalter, R., Myhre, C. E. L., Tsigaridis, K., Vignati, E., Stephanou, E. G., and Wilson, J.: Organic aerosol and global climate modelling: a review, *Atmos. Chem. Phys.*, 5, 1053-1123, 10.5194/acp-5-1053-2005, 2005.

Kim, H., Barkey, B., and Paulson, S. E.: Real refractive indices of α - and β -pinene and toluene secondary organic aerosols generated from ozonolysis and photo-oxidation, *Journal of Geophysical Research: Atmospheres*, 115, <https://doi.org/10.1029/2010JD014549>, 2010.

Krieger, U. K., Siegrist, F., Marcolli, C., Emanuelsson, E. U., Gøbel, F. M., Bilde, M., Marsh, A., Reid, J. P., Huisman, A. J., Riipinen, I., Hyttinen, N., Myllys, N., Kurtén, T., Bannan, T., Percival, C. J., and Topping, D.: A reference data set for validating vapor pressure measurement techniques: homologous series of polyethylene glycols, *Atmos. Meas. Tech.*, 11, 49-63, 10.5194/amt-11-49-2018, 2018.

Kurtén, T., Tiusanen, K., Roldin, P., Rissanen, M., Luy, J.-N., Boy, M., Ehn, M., and Donahue, N.: α -Pinene Autoxidation Products May Not Have Extremely Low Saturation Vapor Pressures Despite High O:C Ratios, *The Journal of Physical Chemistry A*, 120, 2569-2582, 10.1021/acs.jpca.6b02196, 2016.

Laothawornkitkul, J., Taylor, J. E., Paul, N. D., and Hewitt, C. N.: Biogenic volatile organic compounds in the Earth system, *The New phytologist*, 183, 27-51, 10.1111/j.1469-8137.2009.02859.x, 2009.

Lee, B. H., Lopez-Hilfiker, F. D., Mohr, C., Kurtén, T., Worsnop, D. R., and Thornton, J. A.: An Iodide-Adduct High-Resolution Time-of-Flight Chemical-Ionization Mass Spectrometer: Application to Atmospheric Inorganic and Organic Compounds, *Environmental Science & Technology*, 48, 6309-6317, 10.1021/es500362a, 2014.

Lee, C. P., Riva, M., Wang, D., Tomaz, S., Li, D., Perrier, S., Slowik, J. G., Bourgain, F., Schmale, J., Prevot, A. S. H., Baltensperger, U., George, C., and El Haddad, I.: Online Aerosol Chemical Characterization by Extractive Electrospray Ionization-Ultrahigh-Resolution Mass Spectrometry (EESI-Orbitrap), *Environmental Science & Technology*, 54, 3871-3880, 10.1021/acs.est.9b07090, 2020.

Li, H., Almeida, T. G., Luo, Y., Zhao, J., Palm, B. B., Daub, C. D., Huang, W., Mohr, C., Krechmer, J. E., Kurtén, T., and Ehn, M.: Fragmentation inside proton-transfer-reaction-based mass spectrometers limits the detection of ROOR and ROOH peroxides, *Atmos. Meas. Tech.*, 15, 1811-1827, 10.5194/amt-15-1811-2022, 2022.

Li, X., Li, Y., Lawler, M. J., Hao, J., Smith, J. N., and Jiang, J.: Composition of Ultrafine Particles in Urban Beijing: Measurement Using a Thermal Desorption Chemical Ionization Mass Spectrometer, *Environmental Science & Technology*, 55, 2859-2868, 10.1021/acs.est.0c06053, 2021.

Li, Y., Pöschl, U., and Shiraiwa, M.: Molecular corridors and parameterizations of volatility in the chemical evolution of organic aerosols, *Atmos. Chem. Phys.*, 16, 3327-3344, 10.5194/acp-16-3327-2016, 2016.

Lopez-Hilfiker, F. D., Iyer, S., Mohr, C., Lee, B. H., D'Ambro, E. L., Kurtén, T., and Thornton, J. A.: Constraining the sensitivity of iodide adduct chemical ionization mass spectrometry to multifunctional organic molecules using the collision limit and thermodynamic stability of iodide ion adducts, *Atmos. Meas. Tech.*, 9, 1505-1512, 10.5194/amt-9-1505-2016, 2016.

Lopez-Hilfiker, F. D., Pospisilova, V., Huang, W., Kalberer, M., Mohr, C., Stefenelli, G., Thornton, J. A., Baltensperger, U., Prevot, A. S. H., and Slowik, J. G.: An extractive electrospray ionization time-of-flight mass spectrometer (EESI-TOF) for online measurement of atmospheric aerosol particles, *Atmos. Meas. Tech.*, 12, 4867-4886, 10.5194/amt-12-4867-2019, 2019.

Lopez-Hilfiker, F. D., Mohr, C., Ehn, M., Rubach, F., Kleist, E., Wildt, J., Mentel, T. F., Lutz, A., Hallquist, M., Worsnop, D., and Thornton, J. A.: A novel method for online analysis of gas and particle

composition: description and evaluation of a Filter Inlet for Gases and AEROSols (FIGAERO), *Atmos. Meas. Tech.*, 7, 983-1001, 10.5194/amt-7-983-2014, 2014.

Mellouki, A., Wallington, T. J., and Chen, J.: Atmospheric Chemistry of Oxygenated Volatile Organic Compounds: Impacts on Air Quality and Climate, *Chemical Reviews*, 115, 3984-4014, 10.1021/cr500549n, 2015.

Müller, M., Eichler, P., D'Anna, B., Tan, W., and Wisthaler, A.: Direct Sampling and Analysis of Atmospheric Particulate Organic Matter by Proton-Transfer-Reaction Mass Spectrometry, *Analytical Chemistry*, 89, 10889-10897, 10.1021/acs.analchem.7b02582, 2017.

Pai, S. J., Heald, C. L., Pierce, J. R., Farina, S. C., Marais, E. A., Jimenez, J. L., Campuzano-Jost, P., Nault, B. A., Middlebrook, A. M., Coe, H., Shilling, J. E., Bahreini, R., Dingle, J. H., and Vu, K.: An evaluation of global organic aerosol schemes using airborne observations, *Atmos. Chem. Phys.*, 20, 2637-2665, 10.5194/acp-20-2637-2020, 2020.

Peng, Y., Wang, H., Gao, Y., Jing, S., Zhu, S., Huang, D., Hao, P., Lou, S., Cheng, T., Huang, C., and Zhang, X.: Real-time measurement of phase partitioning of organic compounds using a proton-transfer-reaction time-of-flight mass spectrometer coupled to a CHARON inlet, *Atmos. Meas. Tech.*, 16, 15-28, 10.5194/amt-16-15-2023, 2023.

Peräkylä, O., Riva, M., Heikkinen, L., Quéléver, L., Roldin, P., and Ehn, M.: Experimental investigation into the volatilities of highly oxygenated organic molecules (HOMs), *Atmos. Chem. Phys.*, 20, 649-669, 10.5194/acp-20-649-2020, 2020.

Riva, M., Brüggemann, M., Li, D., Perrier, S., George, C., Herrmann, H., and Berndt, T.: Capability of CI-Orbitrap for Gas-Phase Analysis in Atmospheric Chemistry: A Comparison with the CI-API-TOF Technique, *Analytical Chemistry*, 92, 8142-8150, 10.1021/acs.analchem.0c00111, 2020.

Riva, M., Ehn, M., Li, D., Tomaz, S., Bourgain, F., Perrier, S., and George, C.: CI-Orbitrap: An Analytical Instrument To Study Atmospheric Reactive Organic Species, *Analytical Chemistry*, 91, 9419-9423, 10.1021/acs.analchem.9b02093, 2019a.

Riva, M., Rantala, P., Krechmer, J. E., Peräkylä, O., Zhang, Y., Heikkinen, L., Garmash, O., Yan, C., Kulmala, M., Worsnop, D., and Ehn, M.: Evaluating the performance of five different chemical ionization techniques for detecting gaseous oxygenated organic species, *Atmos. Meas. Tech.*, 12, 2403-2421, 10.5194/amt-12-2403-2019, 2019b.

Schobesberger, S., D'Ambro, E. L., Lopez-Hilfiker, F. D., Mohr, C., and Thornton, J. A.: A model framework to retrieve thermodynamic and kinetic properties of organic aerosol from composition-resolved thermal desorption measurements, *Atmos. Chem. Phys.*, 18, 14757-14785, 10.5194/acp-18-14757-2018, 2018.

Shilling, J. E., Chen, Q., King, S. M., Rosenoern, T., Kroll, J. H., Worsnop, D. R., DeCarlo, P. F., Aiken, A. C., Sueper, D., Jimenez, J. L., and Martin, S. T.: Loading-dependent elemental composition of α -pinene SOA particles, *Atmos. Chem. Phys.*, 9, 771-782, 10.5194/acp-9-771-2009, 2009.

Smith, J. N., Moore, K. F., McMurtry, P. H., and Eisele, F. L.: Atmospheric Measurements of Sub-20 nm Diameter Particle Chemical Composition by Thermal Desorption Chemical Ionization Mass Spectrometry, *Aerosol Science and Technology*, 38, 100-110, 10.1080/02786820490249036, 2004.

Song, M., He, S., Li, X., Liu, Y., Lou, S., Lu, S., Zeng, L., and Zhang, Y.: Optimizing the iodide-adduct chemical ionization mass spectrometry (CIMS) quantitative method for toluene oxidation intermediates: experimental insights into functional-group differences, *Atmos. Meas. Tech.*, 17, 5113-5127, 10.5194/amt-17-5113-2024, 2024.

Stark, H., Yatavelli, R. L. N., Thompson, S. L., Kang, H., Krechmer, J. E., Kimmel, J. R., Palm, B. B., Hu, W., Hayes, P. L., Day, D. A., Campuzano-Jost, P., Canagaratna, M. R., Jayne, J. T., Worsnop, D. R., and Jimenez, J. L.: Impact of Thermal Decomposition on Thermal Desorption Instruments: Advantage of Thermogram Analysis for Quantifying Volatility Distributions of Organic Species, *Environmental Science & Technology*, 51, 8491-8500, 10.1021/acs.est.7b00160, 2017.

Stein, S. E. and Scott, D. R.: Optimization and testing of mass spectral library search algorithms for compound identification, *Journal of the American Society for Mass Spectrometry*, 5, 859-866, 10.1016/1044-0305(94)87009-8, 1994.

Thornton, J. A., Mohr, C., Schobesberger, S., D'Ambro, E. L., Lee, B. H., and Lopez-Hilfiker, F. D.: Evaluating Organic Aerosol Sources and Evolution with a Combined Molecular Composition and Volatility Framework Using the Filter Inlet for Gases and Aerosols (FIGAERO), *Accounts of Chemical Research*, 53, 1415-1426, 10.1021/acs.accounts.0c00259, 2020.

Tikkanen, O. P., Buchholz, A., Ylisirniö, A., Schobesberger, S., Virtanen, A., and Yli-Juuti, T.: Comparing secondary organic aerosol (SOA) volatility distributions derived from isothermal SOA particle evaporation data and FIGAERO-CIMS measurements, *Atmos. Chem. Phys.*, 20, 10441-10458, 10.5194/acp-20-10441-2020, 2020.

Vogel, A. L., Schneider, J., Müller-Tautges, C., Klimach, T., and Hoffmann, T.: Aerosol Chemistry Resolved by Mass Spectrometry: Insights into Particle Growth after Ambient New Particle Formation, *Environmental Science & Technology*, 50, 10814-10822, 10.1021/acs.est.6b01673, 2016.

Voisin, D., Smith, J. N., Sakurai, H., McMurry, P. H., and Eisele, F. L.: Thermal Desorption Chemical Ionization Mass Spectrometer for Ultrafine Particle Chemical Composition, *Aerosol Science and Technology*, 37, 471-475, 10.1080/027868203000959, 2003.

Wang, D. S., Lee, C. P., Krechmer, J. E., Majluf, F., Tong, Y., Canagaratna, M. R., Schmale, J., Prévôt, A. S. H., Baltensperger, U., Dommen, J., El Haddad, I., Slowik, J. G., and Bell, D. M.: Constraining the response factors of an extractive electrospray ionization mass spectrometer for near-molecular aerosol speciation, *Atmos. Meas. Tech.*, 14, 6955-6972, 10.5194/amt-14-6955-2021, 2021.

Wilson, J., Imre, D., Beránek, J., Shrivastava, M., and Zelenyuk, A.: Evaporation Kinetics of Laboratory-Generated Secondary Organic Aerosols at Elevated Relative Humidity, *Environmental Science & Technology*, 49, 243-249, 10.1021/es505331d, 2015.

Witkowski, B., al-Sharafi, M., Błaziak, K., and Gierczak, T.: Aging of α -Pinene Secondary Organic Aerosol by Hydroxyl Radicals in the Aqueous Phase: Kinetics and Products, *Environmental Science & Technology*, 57, 6040-6051, 10.1021/acs.est.2c07630, 2023.

Wong, C., Liu, S., and Nizkorodov, S. A.: Highly Acidic Conditions Drastically Alter the Chemical Composition and Absorption Coefficient of α -Pinene Secondary Organic Aerosol, *ACS Earth and Space Chemistry*, 6, 2983-2994, 10.1021/acsearthspacechem.2c00249, 2022.

Xu, L., Coggon, M. M., Stockwell, C. E., Gilman, J. B., Robinson, M. A., Breitenlechner, M., Lamplugh, A., Crounse, J. D., Wennberg, P. O., Neuman, J. A., Novak, G. A., Veres, P. R., Brown, S. S., and Warneke, C.: Chemical ionization mass spectrometry utilizing ammonium ions (NH_4^+ CIMS) for measurements of organic compounds in the atmosphere, *Atmos. Meas. Tech.*, 15, 7353-7373, 10.5194/amt-15-7353-2022, 2022.

Yang, L. H., Masayuki, T., Yunle, C., and Ng, N. L.: Characterization of thermal decomposition of oxygenated organic compounds in FIGAERO-CIMS, *Aerosol Science and Technology*, 55, 1321-1342, 10.1080/02786826.2021.1945529, 2021.

Ylisirniö, A., Barreira, L. M. F., Pullinen, I., Buchholz, A., Jayne, J., Krechmer, J. E., Worsnop, D. R., Virtanen, A., and Schobesberger, S.: On the calibration of FIGAERO-ToF-CIMS: importance and impact of calibrant delivery for the particle-phase calibration, *Atmos. Meas. Tech.*, 14, 355-367, 10.5194/amt-14-355-2021, 2021.

Zaytsev, A., Breitenlechner, M., Koss, A. R., Lim, C. Y., Rowe, J. C., Kroll, J. H., and Keutsch, F. N.: Using collision-induced dissociation to constrain sensitivity of ammonia chemical ionization mass spectrometry (NH_4^+ CIMS) to oxygenated volatile organic compounds, *Atmos. Meas. Tech.*, 12, 1861-1870, 10.5194/amt-12-1861-2019, 2019.

Zhang, X., McVay, R. C., Huang, D. D., Dalleska, N. F., Aumont, B., Flagan, R. C., and Seinfeld, J. H.: Formation and evolution of molecular products in α -pinene secondary organic aerosol, *Proceedings of the National Academy of Sciences*, 112, 14168-14173, doi:10.1073/pnas.1517742112, 2015.

Zhao, J., Häkkinen, E., Graeffe, F., Krechmer, J. E., Canagaratna, M. R., Worsnop, D. R., Kangasluoma, J., and Ehn, M.: A combined gas- and particle-phase analysis of highly oxygenated organic molecules (HOMs) from α -pinene ozonolysis, *Atmos. Chem. Phys.*, 23, 3707-3730, 10.5194/acp-23-3707-2023, 2023.

803 Zhao, J., Mickwitz, V., Zhang, J., Alton, M., Canagaratna, M., Graeffe, F., Schobesberger, S., Worsnop,
804 D., and Ehn, M.: Comparison of Gaseous and Particulate Highly Oxygenated Organic Molecules from
805 the Ozonolysis of Terpenes, *ACS ES&T Air*, 1, 1294-1303, 10.1021/acsestair.4c00121, 2024a.
806 Zhao, J., Mickwitz, V., Luo, Y., Häkkinen, E., Graeffe, F., Zhang, J., Timonen, H., Canagaratna, M.,
807 Krechmer, J. E., Zhang, Q., Kulmala, M., Kangasluoma, J., Worsnop, D., and Ehn, M.: Characterization
808 of the Vaporization Inlet for Aerosols (VIA) for online measurements of particulate highly oxygenated
809 organic molecules (HOMs), *Atmos. Meas. Tech.*, 17, 1527-1543, 10.5194/amt-17-1527-2024, 2024b.
810 Zuth, C., Vogel, A. L., Ockenfeld, S., Huesmann, R., and Hoffmann, T.: Ultrahigh-Resolution Mass
811 Spectrometry in Real Time: Atmospheric Pressure Chemical Ionization Orbitrap Mass Spectrometry of
812 Atmospheric Organic Aerosol, *Analytical Chemistry*, 90, 8816-8823, 10.1021/acs.analchem.8b00671,
813 2018.

814